Invacare® Standing Transfer Vest





EN	Sling	User Manual	1
DA	Sejl	Brugsanvisning	3
DE	Patientengurte	Gebrauchsanweisung	4
ES	Eslinga	Manual del usuario	5
FR	Sangles	Manuel d'utilisation	6
IT	Imbracatura	Manuale d'uso	8
NL	Liftmat	Gebruiksaanwijzing	9
NO	Seil	Bruksanvisning	10
PT	Cesta de transferência	Manual de utilização	11
SV	Lyftselar	Bruksanvisning	12

EN

I General

I.1 Figures

Figur 1
 Figur
 Figur
 Abbildung
 Afbeelding
 Figures
 Figura
 Figura
 Figura

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

1.2 Symbols

	Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words. WARNING! Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury. CAUTION! Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.
	Gives useful tips, recommendations and information for efficient, trouble-free use.
	IMPORTANT indicates a hazardous situation that could result in damage to property if it is not avoided.
	This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.
	Safe Working Load
	Consult the user manual.
	DO NOT exceed 95°C washing temperature.
	DO NOT bleach.
	Tumble Dry, Normal, Low Heat
	DO NOT iron.

1.3 Intended Use

The Invacare® Standing Transfer Vest is available with or without groin band, the groin band can be attached to the transfer standing vest very easily.

The Standing Transfer Vest is suitable for clients with reduced standing/walking abilities, and/or reduced balance. The sling provides maximum personal access for dressing and hygiene.

-  The Standing Transfer Vest is applied around the upper part of the body, located in the under arm area. It is essential that the user has the ability to stand up, and has some degree of weight bearing ability to ensure that he feels contact with the floor.

1.4 Service Life

All Invacare slings are guaranteed for one year. Expected service life of the product is 1-5 years. Service life time varies depending on the fabric, frequency of use, laundering practices, and weight carried.

2 Safety

2.1 General Guidelines

WARNING!

- DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.
- Please check your sling on delivery. If the sling is not in perfect condition, please contact your supplier.

2.2 Safe Working Load

	Solid	
XS	1485404	200 kg
S	1485405	200 kg
M	1485406	200 kg
L	1484683	200 kg
XL	1484684	200 kg

3 Usage

3.1 Sling application

WARNING!

- Slings should be checked regularly for wear and damage (in accordance with each country's regulations). NEVER use a damaged sling.
- If you are in any doubt, please do not attempt to use the sling and seek the advice of a professional or your supplier.
- Invacare recommends that a risk assessment is always carried out by a professional prior to issuing lifting equipment. Always choose the sling design and size according to the client's weight, size and physical ability.
- If you are in doubt about the use of Invacare's slings with other manufacturers' equipment, contact Invacare or your supplier for advice.

Sling application — Sitting in chair or wheelchair

1. Lean the client forward ensuring that he is well supported to reduce risk of falling and to promote feelings of security.
2. Slide the vest down behind the back of the client in a central position aligned with the waist.
3. Sit the client back against the sling. **Figure 1.**
4. Take the arms through the straps and fasten the vest around the client's stomach by using the adjustable buckle. **Figure 2.**

3.2 Attach loops to the spreader bar or lifting strap

WARNING!

- Use only CE marked spreader bars and hoists, preferably with a rated Safe Working Load equal to or higher than that of the sling. If the SWL is lower than of the sling, the SWL of the spreader bar or hoist is applicable.
- Always check that all the straps are properly attached to the spreader bar or lifting straps, before attempting to lift the client.

1. Attach the loops to the spreader bar or lifting straps of the lifting equipment.

The sling can be used with spreader bars with 2 attachment points.

For security and comfort reasons, the size (width) of the spreader bar must correspond with the size of the sling.

2. Before lifting the client, always check that both feet are placed on the floor for weight bearing.
3. Help the client to lean forward (in his natural way of movement) before beginning the lift. **Figure 3.**
4. When the client is in a standing position, the posture of the client can be adjusted by altering the shoulder straps. **Figure 4.**

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© 2014 Invacare® Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.
Invacare Corporation