Invacare® InvaSoft® Hybrid

| | | |
|----|---|-----|
| en | Mattress system User Manual | 3 |
| fr | Système de matelas Manuel d'utilisation | 13 |
| da | Madrassystem Brugsanvisning | 25 |
| de | Matratzensystem Gebrauchsanweisung | 35 |
| es | Sistema de colchón Manual del usuario | 45 |
| pt | Sistema de colchões Manual de utilização | 55 |
| it | Sistema materasso Manuale d'uso | 67 |
| nl | Matrassysteem Gebruiksaanwijzing | 77 |
| fi | Patjajärjestelmä Käyttöohje | 87 |
| sv | Madrassystem Bruksanvisning | 97 |
| no | Madrassystem Bruksanvisning | 107 |

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved
for future reference.

Yes, you can.®

©2025 Invacare International GmbH

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare International GmbH or its subsidiaries unless otherwise noted.

Contents

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved
for future reference.

| | |
|--|-----------|
| 1 General | 4 |
| 1.1 Introduction | 4 |
| 1.2 Symbols in this Document | 4 |
| 1.3 Compliance | 4 |
| 1.4 Warranty Information | 4 |
| 1.5 Limitation of Liability | 4 |
| 1.6 Service Life | 4 |
| 2 Safety | 5 |
| 2.1 Safety Information | 5 |
| 2.2 Safety Information for Transport | 5 |
| 2.3 Fire resistance Safety Information | 5 |
| 3 Product Overview | 6 |
| 3.1 Product Description | 6 |
| 3.2 Intended Use | 6 |
| 3.2.1 Indications | 6 |
| 3.2.2 Contraindications | 6 |
| 3.3 Components | 6 |
| 3.4 Symbols on Product | 6 |
| 4 Usage | 8 |
| 4.1 Safety Information | 8 |
| 4.2 Preparing Mattress for Use | 8 |
| 5 Maintenance | 9 |
| 5.1 Inspection | 9 |
| 5.2 Cleaning and Disinfection | 9 |
| 5.2.1 General Safety Information | 9 |
| 5.2.2 Cleaning Intervals | 9 |
| 5.2.3 Cleaning Instructions | 9 |
| 5.2.4 Disinfection Instructions | 9 |
| 5.3 Replacing Cover | 10 |
| 6 After Use | 11 |
| 6.1 Storage | 11 |
| 6.2 Reconditioning | 11 |
| 6.3 Disposal | 11 |
| 7 Technical Data | 12 |
| 7.1 Mattress Specifications | 12 |
| 7.2 Materials | 12 |
| 7.3 Environmental Conditions | 12 |

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Only use this product if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.2 Symbols in this Document

Symbols and signal words are used in this document and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

NOTICE

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

Other Symbols

UK Responsible Person

Indicates if a product is not manufactured in the UK.

Triman

Indicates recycling and sorting rules (only relevant for France).

1.3 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

1.4 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

1.5 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

1.6 Service Life

The expected service life of this product is five years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

2 Safety

2.1 Safety Information

Education, clinical judgement and action-based planning based on vulnerability are fundamental factors in prevention of pressure injuries.

A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

WARNING!

Risk of Serious Injury or Damage

Improper use of this product may cause injury or damage.

- If you are unable to understand the warnings, cautions or instructions, contact a health care professional or provider before attempting to use this equipment.
- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manual or instruction sheet supplied with this product or optional equipment.

WARNING!

Risk of fire or explosion!

A cigarette can burn a hole in the bed surface and cause damage to the device. Also, patient clothing, bed sheets, etc. may be combustible and cause a fire. Failure to observe this warning can result in a severe fire, property damage and cause physical injury or death.

- Take special care in oxygen rich environments.
- Do not smoke.

WARNING!

Risk of Developing Pressure Injuries

The device delivers effective pressure redistribution, when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these is the only item deployed between the support surface and the user.

- Bed sheets must be loosely fitted, with creases smoothed out.
- Ensure that the support surface in contact with the user is kept free from crumbs and other food debris, and that drip cables, stents, and other foreign objects do not become entrapped between the support surface and the user.
- Heated over blankets must only be used in consultation with a qualified healthcare professional, as an increase in temperature can increase the risk of developing pressure injuries.

CAUTION!

Risk of Injury

Non-original or wrong parts may affect the function and safety of the product.

- Only use original parts for the product in use.
- Due to regional differences, refer to your local Invacare catalogue or website for available options or contact your Invacare distributor. See addresses at the end of this document.

2.2 Safety Information for Transport

- Take care when handling the product to ensure no damage to the cover.
- It is recommended that two people lift/carry the product.
- Avoid contact with jewellery, nails, abrasive surfaces, etc.
- Do not drag the product.
- Avoid contact with wall, door frames, door catches or locks, etc.
- Do not transport in roll cages unless completely protected from the sharp edges of the cage.

2.3 Fire resistance Safety Information

CAUTION!

Risk of diminished fire resistance

If you do not use the **Stockinette** as recommended the fire resistance properties of the product are compromised.

- ALWAYS use the **Stockinette** to encapsulate all mattress components within the cover.

3 Product Overview

3.1 Product Description

The Invacare InvaSoft® Hybrid mattress acts as a reactive pressure redistributing support/mattress for patients at High Risk that can, by the airflow throughout the air cells, promote a pressure distribution support surface.

Due to its unique construction combining foam, noodle material, and air, this mattress provides high user comfort, as well as a solution for pressure care and microclimate management without the need for an external pump or electricity. The system uses the user weight and its movements to distribute the internal air reservoir amongst the cells, to always maintain optimal pressure redistribution. The mattress has a sloped heel section for protection of the vulnerable heels.

The high-breathability water-resistant cover provides a vapor-permeable, multi-stretch surface, with enhanced microclimate properties, to promote patient comfort, contribute to skin integrity and maximize the effectiveness of the underlying layers.

The InvaSoft® mattress is produced in high-quality materials for function and durability.

3.2 Intended Use

This pressure redistribution mattress is intended to be used in conjunction with an appropriately sized bed frame and side rails. Check the bed's user manual for the appropriate size of the mattress to be used.

This product delivers effective pressure redistribution to users, when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these is the only item deployed between the support surface and the user.

Intended Users

Adults and adolescents limited in their ability to change and control body position, and/or with limited or no physical activity in a hospital, long term care or home care setting.

3.2.1 Indications

Suitable for supporting the management of all categories of pressure injuries when combined with an individual and comprehensive pressure injury protocol.

Suitable for use in home care, residential, nursing and acute care settings.

3.2.2 Contraindications

Not suitable for users with unstable spinal cord fractures and/or cervical traction.

Always consult a physician prior to using this device.

3.3 Components

The following components are included within the scope of delivery:


Fig. 3-1

- Ⓐ Top cover
 - Ⓑ Non-woven top layer
 - Ⓒ Noodle layer
 - Ⓓ Headrest foam
 - Ⓔ Internal spacers
 - Ⓕ Air cell assembly
 - Ⓖ Anti-shear material
 - Ⓗ Heel sloped foam
 - Ⓘ U-core foam
 - ⓵ Bottom cover
 - Ⓚ User Manual
- Stockinette encapsulating all mattress components within the cover (not shown)

3.4 Symbols on Product

| | | | |
|--|-------------------------|--|---|
| | European Conformity | | UK Conformity Assessed |
| | European representative | | Medical Device |
| | Manufacturer | | Date of manufacture |
| | Read user manual | | User weight limit. See 7 Technical Data, page 12. |
| | LOT number | | Do not put near flame |

| | | | |
|---|--|---|---------------|
|  | Do not pierce or cut |  | Hand wash |
|  | Machine wash temperature. Max. temperature is shown in symbol. |  | Do not bleach |

| | | | |
|--|------------------|---|----------------------------|
|  | Do not dry clean |  | Tumble dry low temperature |
|  | Line dry |  | Do not iron |
|  | | Orientation symbol | |

4 Usage

4.1 Safety Information

WARNING!

Insufficient repositioning of the patient can lead to tissue compression and potential injury formation.

To relieve pressure, it is very important for the patients to reposition themselves, or to be repositioned, on a regular basis.

- Before using the product, always consult a qualified healthcare professional for clinical judgement.
- Monitor the patient frequently.

CAUTION!

- Make sure that the printed side of the mattress cover always faces upwards.
- Ensure that the distance between the surface of the mattress and the top of the side rail is at least 220 mm.

CAUTION!

Risk of damage to the mattress

If holes are present in the mattress cover, there is a risk that liquids may ingress and contamination may occur.

- Ensure that the mattress is not jammed or damaged by sharp edges.
- Do not place hypodermic needles, venflons, scalpels or other similarly sharp objects on or under the mattress.
- Do not use electrically heated blankets directly on or under the mattress.
- Ensure that all venflons are taped down correctly with no sharp edges exposed.
- When using bridging boards or other patient transfer aids, check for any sharp edges or burrs before use.
- When using the mattress on a profiling bed ensure that the knee break is used before the backrest.
- Attach medical equipment including infusion pumps and monitors to appropriate bed accessories.
- Avoid cigarette burns and pet claws puncturing the mattress cover.

4.2 Preparing Mattress for Use

CAUTION!

Risk of Injury or Damage to Property

- Before use, check the product for any defect.
- Do not use a defective product.

1. Remove all packaging before use.
2. Place mattress directly on bed frame.
3. Make sure that heel slope section is oriented at foot end of bed frame.

5 Maintenance

5.1 Inspection

It is recommended to check mattresses (air cells and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the release of each patient or after each period of use by a suitably qualified and competent person.

Check mattresses

1. Unzip the cover completely.
2. Check for any staining on the white underside of the cover.
3. Check for any staining on the internal elements.
4. Replace any stained items and dispose of as per local authority procedure.

5.2 Cleaning and Disinfection

5.2.1 General Safety Information

CAUTION!

Risk of Contamination

- Take precautions for yourself and use appropriate protective equipment.

NOTICE!

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking it into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

5.2.2 Cleaning Intervals

NOTICE!

Regular cleaning and disinfection enhances smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product

- regularly while it is in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

NOTICE!

- Keep a cleaning record as part of cleaning the system.

5.2.3 Cleaning Instructions

NOTICE!

- The product does not tolerate cleaning in automatic washing plants, with high-pressure cleaning equipment or steam.

A cleaning record should be kept as part of cleaning the system.

Cleaning covers

(Removal of contaminants such as dust and organic matter)

1. Remove all covers for laundering.
2. Launder covers with maximum temperature as stated on product label, using a diluted detergent solution (Instructions on label).

NOTICE!

- Washing at higher temperatures will cause shrinkage.

Drying covers

1. Hang covers from a line or bar and drip dry in a clean indoor environment or tumble dry on a low heat setting.

NOTICE!

- Tumble dry setting must not exceed 40 °C.
- Do not tumble dry for longer than 10 minutes.
- Dry thoroughly before use.

5.2.4 Disinfection Instructions

NOTICE!

- Only use disinfectants and methods approved by your local infection control institution and follow your local infection control policy.
- Follow your local decontamination protocols.

Disinfecting Covers

(Reducing the number of microorganisms)

NOTICE!

- Failure in disinfection process may result in the accumulation of reagent that could damage the polyurethane coating, react with the bed frame, or negate the biocompatibility results.
- Ensure that all cleaning agents, and disinfectants, are thoroughly rinsed off.
- Dry thoroughly before use.

1. Wipe down the cover with a suitable detergent.
2. Rinse the cover thoroughly with clean water using a single use nonabrasive cloth.
3. Dry the cover thoroughly.

NOTICE!

Heavy soilage

Where the mattress is badly soiled, we recommend cleaning with a dilute cleaning solution at maximum temperature in the washing machine (see product label).

- Clean up all spillages of bodily fluids i.e. blood, urine, faeces, sputum, wound exudate and all other bodily secretions as soon as possible using a suitable detergent.
- Large spillages of blood should be absorbed and removed with paper towels first.

WARNING!

- In event of contamination, contact your hygiene specialist.
- Remove contaminated foams from use.

Autoclaving Covers

Autoclave the cover at 110 °C.

NOTICE!

- Make sure the cover is outstretched (not folded).
- Do not place covers on top of each other.

5.3 Replacing Cover

1. Unzip cover and remove it carefully from foam insert.
2. Place new cover onto foam insert.
3. Close zipper.

NOTICE!

- Ensure that the corners of the foam insert are positioned correctly into the corners of the cover.
- Ensure that the castellated foam is facing upwards when inserted into its cover.

WARNING!

- In event of contamination, contact your hygiene specialist.
- Remove contaminated foams from use.

6 After Use

6.1 Storage

NOTICE!

- Store mattresses in a dry environment.
- Store mattresses within a protective cover.
- Store items on clean, dry (off-floor) surface with no sharp edges to avoid any possible damage.
- Never store other items on top of a mattress.
- Do not store mattresses next to radiators or other heating devices.
- Protect mattresses from direct sunlight.

For environmental conditions for storage, see chapter 7 *Technical Data*, page 12.

6.2 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection
- Cleaning and disinfection

For detailed information, see 5 *Maintenance*, page 9.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

6.3 Disposal

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

7 Technical Data

7.1 Mattress Specifications

| Product | Dimensions [mm] | | | Air Cell Height [mm] | Maximum User Weight [kg] | Weight of Product [kg] ¹ |
|-----------------|-----------------|------------------|--------|----------------------|--------------------------|-------------------------------------|
| | Length | Width | Height | | | |
| InvaSoft Hybrid | 830, 850, 880 | 1970, 2000, 2100 | 170 | 50 | 230 | 18 |

1 Based on the weight of a standard size mattress. This can change if different sizes are ordered.

7.2 Materials

| | |
|--|--|
| Stockinette | 65% modified acrylic, 35% cotton |
| Cover | Polyurethane coated nylon knitted fabric |
| Foam (U-core, headrest, heel slope) | Polyurethane foam |
| Noodle layer, internal spacer | Polyester |
| Top layer | Carbide fiber non woven |

All product components do not contain natural rubber latex.

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Invacare AB
Invacare Australia Pty. Ltd.
Invacare UK Operations Limited