

Invacare® | *Softform® Heelpad*



Pressure Reducing Heelpad

The *Invacare Softform Heelpad* is designed to reduce pressure on the vulnerable heel area. The controlled volume of gel in the sacs allows pressure displacement and even weight distribution. Pressure on a patient's heels are kept to a minimum when the heels are in contact with the gel sacs. The *Invacare Softform Heelpad* was designed to be used in conjunction with the *Invacare Odstock Wedge*, but can also be used independently as an effective pressure reducing device.

Softform®



Yes, you can.®



Features and Options

Designed to reduce pressure on the heels

The *Invacare Softform Heelpad* enables vulnerable heels to be immersed into independent gel sacs, offering significant pressure reduction and comfort to clients.

Enhanced combined performance

The *Softform Heelpad* offers enhanced performance, when used in conjunction with the *Softform Odstock Wedge*, by providing improved support for the lower limbs, and pressure reduction for the heels.



Single and Double Heelpads available

Single and Double formats available to meet specific patient requirements.



Heelpads supplied with three platilon covers

Platilon protective covers are disposable, and help to avoid contamination of the heelpad.



Designed to be used in conjunction with the Invacare Odstock Wedge

Used in conjunction, these two products combine to significantly reduce pressure on the vulnerable heel area.

Technical data

Softform Heelpad					
	Width: Double: 620mm Single: 320mm	Length: 480 mm	Height: 50 mm, 35 mm	Weight: 3.58 kg	Cleaning: 80°

Manufactured to European Standards:
Interior Foam: BS 5852: 1982 Crib 5, BSEN 597-1,2
Cover: BS 7175: 1989 Crib 5

Colour



Fuchsia

Warranty 1 Year

